

Best-in-class optimized biopreservation media for cells and tissues

- 27% DMSO in Saline
- Sterile Contents
- 80mL Fill Volume
- USP/Highest Quality Components
- cGMP Manufactured
- Sterility, Endotoxin, and DMSO Content Testing



BloodStor®, is a series of generic cGMP freeze media products used to cryopreserve stem cells, other cells, and cell components isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor® products are formulated using USP grade components and Water for Injection (WFI) quality water.

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

Glossary of label symbols:

REF	LOT	STERILE
Order Number	Number Batch Code	Aseptic Processing Technique
Exp. Date	Temperature	Manufactured by



3303 Monte Villa Parkway, Suite 310
Bothell, WA 98021 USA
1.866.424.6543 Phone | 1.425.402.1433
BioLifeSolutions.com

Ordering Information

Product Name	Size	Part #
BloodStor® 27 NaCl	80mL bag	327207

Related Products

BloodStor® 55-5, BloodStor® 100
Please visit website for ordering information.

Shipping and Storage

BloodStor® 27 NaCl ships at ambient temperature.
Upon receipt, store at 2° to 25° C, protected from light,
until ready to use.

To Order

Call: 1.866.424.6543 | **Fax:** 425.402.1433
Sales: SalesOne@BioLifeSolutions.com
Web: BioLifeSolutions.com
Technical Support: info@BioLifeSolutions.com