

## Best-in-class optimized biopreservation media for cells and tissues

- 55% (w/v) DMSO
- 5% (w/v) Dextran 40
- Standard and Custom Fill Volumes
- USP/Highest Quality Components
- cGMP Manufactured
- Sterility, Endotoxin, and Cell-Based Release Testing



BloodStor<sup>®</sup> is a series of generic cGMP freeze media products used to cryopreserve stem cells, other cells, and cell components isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor<sup>®</sup> products are formulated using USP grade components and Water for Injection (WFI) quality water.

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

#### Glossary of label symbols:

 REF	 LOT	 STERILE
Order Number	Number Batch Code	Aseptic Processing Technique
 Exp. Date	 Temperature	 Manufactured by



3303 Monte Villa Parkway, Suite 310  
Bothell, WA 98021 USA  
1.866.424.6543 Phone | 1.425.402.1433  
BioLifeSolutions.com

### Ordering Information

Product Name	Size	Part #
BloodStor <sup>®</sup> 55-5	5.6mL vial	355371
BloodStor <sup>®</sup> 55-5	7.2mL vial	355372

### Related Products

BloodStor<sup>®</sup> 100  
Please visit website for ordering information.

### Shipping and Storage

BloodStor<sup>®</sup> 55-5 ships at ambient temperature.  
Upon receipt, store at 2° to 8° C, until ready to use.

### To Order

**Call:** 1.866.424.6543 | **Fax:** 425.402.1433  
**Sales:** SalesOne@BioLifeSolutions.com  
**Web:** BioLifeSolutions.com  
**Technical Support:** info@BioLifeSolutions.com