

CRYOPRESERVE CELLS

BloodStor® 27 NaCl

Preserving the promise™


Best-in-class optimized biopreservation media for cells and tissues

- ◆ 27% DMSO in Saline
- ◆ Sterile Contents
- ◆ USP/Highest Quality Components
- ◆ cGMP Manufactured
- ◆ Sterility, Endotoxin, and DMSO Content Testing

BloodStor®, is a series of generic cGMP freeze media products used to cryopreserve stem cells, other cells, and cell components isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. All BloodStor® products are formulated using USP grade components and Water for Injection (WFI) quality water.

Test methods and criteria are provided on all lot specific Certificates of Analysis and Release.

Glossary of label symbols

REF Order Number	LOT Number Batch Code	STERILE Aseptic Processing Technique
 Exp. Date	 Temperature	 Manufactured by

Ordering Information

Product Name	Size	Part #
BloodStor® 27 NaCl	80mL bag	327207

Related Products

BloodStor® 55-5, BloodStor® 100
Please visit our website for ordering information.

Shipping and Storage

BloodStor® 27 NaCl ships at ambient temperature.
Upon receipt, store at 2° to 25°C, protected from light, until ready to use.

To Order

Call: 1.866.424.6543 | **Fax:** 1.425.402.1433
Sales: SalesOne@BioLifeSolutions.com
Web: BioLifeSolutions.com
Technical Support: info@BioLifeSolutions.com